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4. Holden LK, Finley CC, Firszt JB, Holden TA, Brenner C, Potts LG, et al. Factors affecting open-set word recognition in adults with cochlear implants. Ear Hear.2013 May-Jun;34(3):342-60.

5. Stakhovskaya , Sridhar D, Bonham BH, Leake PA. Frequency map for the human cochlear spiral ganglion: implications for cochlear implants. J Assoc Res Otolaryngol. 2007 jun;8(2):220-33.

6. Ariyasu L, Galey FR, Hilsinger R, Jr., Byl FM. Computer-generated three-dimensional reconstruction of the cochlea. Otolaryngol Head Neck Surg. 1989 Feb;100(2):87-91.

7. Roland JT, Jr. A model for cochlear implant electrode insertion and force evaluation: results with a new electrode design and insertion technique. Laryngoscope. 2005 Aug;115(8):1325-39

8. Aschendorff A, Kromeier J, Klenzner T, Laszig R. Quality control after insertion of the Nucleus@ Contour@ and Contour Advance@electrode in adults. Ear Hear. 2007 Apr;28(2 Suppl):75S-9S.

9. Balkany T, Hodges A, Menapace C, Hazard L, Driscoll C, Gantz B, et al. Nucleus Freedom North American clinical trial. Otolaryngol Head Neck Surg. 2007 May;136(5):757-62.

10. Mukherjee P, Uzun-Coruhlu H, Wong CC, Curthoys IS, Jones AS, Gibson WP. Assessment of intracochlear trauma caused by the insertion of a new straight research array. Cochlear Implants Int. 2012 Aug;13(3):156-62.

11. Skarzynski H, Lorens A, Matusiak M, Porowski M, Skarzynski PH, James CJ. Cochlear implantation with the Nucleus@ Slim Straight electrode in subjects with residual low-frequency hearing. Ear Hear. 2014 Mar-Apr;35(2):e33-43.

12. Jurawitz MC, Buchner A, Harpel T, Schussler M, Majdani O, Lesinski-Schiedat A, et al. Hearing preservation outcomes with different cochlear implant electrodes: Nucleus@ Hybrid™ L24 and Nucleus Freedom™ C1422. Audiol Neurotol. 2014;19(5):293-309.

13. Aschendorff A, Jaekel K, Klenzner T, Laszig R. Impact of electrode design on facial nerve stimulation in otosclerosis. Cochlear Implants Int. 2004 Sep;5 Suppl 1:63-5.

14. Papsin BC. Cochlear implantation in children with anomalous cochleovestibular anatomy. Laryngoscope. 2005 Jan;115(1 Pt 2 Suppl 106):1-26.

15. Rinia AB, Olphen AF, Dunnebier EA. Cochlear implantation in obstructed cochleas: the effect of the degree of obstruction on the number of activated electrodes and the amount of postoperative speech perception. Clin Otolaryngol. 2006 Aug;31(4):280-6.

16. Coelho DH, Waltzman SB, Roland JT, Jr. Implanting common cavity malformations using intraoperative fluoroscopy. Otol Neurotol. 2008 Oct;29(7):914-9.

17. Tysome JR, Axon PR, Donnelly NP, Evans DG, Ferner RE, O'Connor AF, et al. English consensus protocol evaluating candidacy for auditory brainstem and cochlear implantation in neurofibromatosis type 2. Otol Neurotol. 2013 Dec;34(9):1743-7.

18. Sennaroglu L, Sarac S, Ergin T. Surgical results of cochlear implantation in malformed cochlea. Otol Neurotol. 2006 Aug;27(5):615-23.

19. Sennaroglu L, Ziyal I, Atas A, Sennaroglu C, Yucel E, Sevinc S, et al. Preliminary results of auditory brainstem implantation in prelingually deaf children with inner ear malformations including severe stenosis of the cochlear aperture and aplasia of the cochlear nerve. Otol Neurotol. 2009 Sep;30(6):708-15.

20. Det™an S, Sadeghi-Barzalighi A, Ambett R, Dowell R, Trotter M, Briggs R. Cochlear implants in forty-eight children with cochlear and/or vestibular abnormality. Audiol Neurotol. 2011;16(4):222-32.

21. Loundon N, Rouillon I, Munier N, Marlin S, Roger G, Garabedian EN. Cochlear implantation in children with internal ear malformations. Otol Neurotol. 2005 Jul;26(4):668-73.

22. Sennaroglu L, Colletti V, Manrique M, Laszig R, Offeciers E, Saeed S, et al. Auditory brainstem implantation in children and non-neurofibromatosis type 2 patients: a consensus statement. Otol Neurotol. 2011 Feb;32(2):187-91.

23. Colletti V, Shannon R, Carner M, Veronese S, Colletti L. Outcomes in nontumor adults fitted with the Auditory Brainstem Implant: 10 years' experience.Otol Neurotol. 2009 Aug;30(5):614-8.

24. Roland JT, Jr., Coelho DH, Pantelides H, Waltzman SB. Partial and Double-Array implantation of the ossified cochlea. Otol Neurotol. 2008 Dec; 29 (8):1068-75.

25. Cochlear Limited. D1805413. Cochlear Nucleus Reliability Report, Volume 19, December 2020.

26. Cochlear Implant Reliability MED-EL [Internet]. Medel.com. 2020 [cited 11 December 2020]. Available from: https://www.medel.com/hearing-solutions/cochlear-implants/reliability.

27. Advanced Bionics Reliability Report 2020. 027-N258-02 Rev A. Advanced Bionics AG and affiliates.; 2020

28. Oticon Medical Reliability Report 2020. In accordance with European and Global Consensus on Cochlear Implant Failures and Explantations. 224811UK - version A / 2020.10. Data valid as of 30 June 2020

29. Oticon Medical Reliability Report 2020. Reporting to ANSI/AAM I C186 Standard. 224812US - version A / 2020.11. Data valid as of 30 June 2020

30. Cochlear Nucleus Profile Plus Series 99.97% combined Cumulative Survival Percentage (CSP) within 2 years. CSP includes both device and accident-related issues.

31. The following implants are designed to be compatible for MRI at 1.5 T and 3.0 T with magnet in place (following the instructions in the MRI Guidelines): Nucleus Profile Plus Series: C1612, C1622, C1632. The following implants are designed to be compatible for MRI at 3.0 T with the magnet removed (should your doctor specifically recommend this procedure) and 1.5 T with magnet in place with use of the Cochlear Nucleus Implant Bandage and Splint Kit (following the instructions in the MRI Guidelines): Nucleus Profile Series: C1512, C1522, C1532. Nucleus C124RE: C1422, C124REH, C124RE(CA), and C124RE(ST)

32. Compared to all currently available receiver stimulators available from Cochlear and other cochlear implant manufacturers Based on published device specification information.

33. Cochlear Nucleus Profile Series 99.79% combined Cumulative Survival Percentage (CSP) within 7 years. CSP includes both device and accident-related issues.

34.Cochlear Nucleus C124RE Series 99.00% combined Cumulative Survival Percentage (CSP) within 16 years. CSP includes both device and accident-related issues.

35. Mauger SJ, Warren C, Knight M, Goorevich M, Nel E. Clinical evaluation of the Nucleus 6 cochlear implant system: performance improvements with SmartSound iQ. International Journal Of Audiology. 2014, Aug; 53(8): 564-576 。[科利耳贊助]。

36.Wolfe J,Parkinson A, Schafer E, Gilden J, Rehwinkel K, Gannaway S, et al. Benefit of a Commercially Available Cochlear Implant Processor With Dual-Microphone Beamforming: A Multi-Center Study. Otology And Neurotology. 2012;33(4): 553-560 。科利耳贊助]。

37. Mauger SJ, et al. Clinical outcomes with the Kanso off-the-ear cochlear implant sound processor. Int J Audiol. 2017, Apr; 56(4): 267-276. [Sponsored by Cochlear].

38. Wolfe J, et al. Benefits of Adaptive Signal Processing in a Commercially Available Cochlear Implant Sound Processor. Otol Neurotol.2015 Aug;36(7):1181-90.

39. Cochlear Ltd. D1660797. CP1150 Sound Processor Interim Clinical Investigation Report. 2020; Jan. Data on file.

40. Cochlear Limited. D1376556. CLTD 5709 Clinical Evaluation of Nucleus 7 Cochlear Implant System. 2018, Jan. Data on file.

41. Wolfe J, et al. Improving hearing performance for cochlear implant recipients with use of a digital, wireless, remote-microphone,audio-streaming accessory. J Am Acad Audiol. 2015 Jun;26(6):532-9.

42. Razza S, et al. Evaluation of speech reception threshold in noise in young Cochlear Nucleus' system 6 implant recipients using two different digital remote microphone technologies and a speech enhancement sound processing algorithm. Int Journal Pediat Otorhinolaryngology, 2017; Oct;103:71-75.

43. Jones M, Hearing in Noise Benefit Using the Cochlear Mini Microphone 2t with Nucleus Sound Processors, C12017 Pediatric 15th Symposium on Cochlear Implants in Children, San Francisco, CA. 2017; July.

44. Cochlear Ltd. D1190805 Sound Processor Size Comparison. 2020; March. Data on file.

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47. Cochlear Ltd. D1710313 CP1150 Battery Life Coverage Technical Report. 2020; Mar. Data on file.

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© Cochlear Limited 2022. D1899160-2022-03